

ZD120E Multi-Parameter Monitor

PLAYBACK CAPABILITIES

- Trend Chart Review: 168h
- Blood Pressure Data: 1500 sets
- Alarm Data Playback: 100 piece
- ECG Waveform Playback: 120min

Simple and Elegant



Easy to use

Designed for the medical staff, easy to move around between wards

High Performance Battery



Only 8 hours to recharge

2 Hours battery operating time

Large Screen



31cm Screen

Color TFT LCD

- Accurate measurement during low perfusion.
- Anti-movement, stable performance.
- Suitable for adult, pediatric and neonatal patients.
- Making use of Motion-tolerant technology - Motion tolerant pulse oximetry is critical for ensuring patient safety because patient movement can thwart accurate readings and delay diagnosis of serious respiratory compromise.
- Provides easy access to critical information.
- Numbers and trend data are provided and clearly displayed in real-time.
- Field-upgradable via a USB port for feature enhancements/improvements or new parameters.

Based on a portable, lightweight design that enables easy transport

Parameters



7 - LEAD ECG



HEART RATE



PULSE RATE




BLOOD PRESSURE



BLOOD OXYGEN



RESPIRATION



BODY TEMPERATURE
(OPTIONAL)


ZD120E MULTI-PARAMETER MONITOR



01

Fanless design

Quiet, energy-saving, clean and quiet, does not affect the patient's rest, and reduces the possibility of cross-infection



02

Sound and light alarm

Three-color warning lights clearly indicate, real-time reminder at any time to prevent dangerous situations from happening



03

Flexible assembly

With different mobile mounting solutions, in line with different environmental needs



04

Expansion interface

Can be connected to different devices to meet different functional requirements of the monitor



REG Number 2021/385832/07 | SAHPRA License: 000001989MD
Corner Barbara and North Reef Road; Elandsfontein, Henville, Germiston,
1429, South Africa

+27 (10) 013 3010 | sales@medhealthsup.com